

Phospho-Alt-Rpl36-S142 Rabbit Polyclonal Antibody - Anti-Alt-Rpl36 S142

Item Number: CM0000566

Introduction

Phospho-specific rabbit polyclonal antibody targeting human alt-RPL36 at serine 142. Validated for Western blot, dot blot, and ELISA applications, suitable for studying phosphorylation-dependent regulatory processes.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-alt-RPL36-S142 Rabbit pAb
Species	Human
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S142 of alt-RPL36
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human
Shipping	Ice bag