

Phospho-Alt-Rpl36-S19 Rabbit Pab - Alt-Rpl36

Item Number: CM0000634

Introduction

Phospho-alt-RPL36-S19 Rabbit Polyclonal Antibody for WB, DB, ELISA. Detects phosphorylated alt-RPL36 at Ser19 in human samples. Suitable for phosphorylation studies and signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-alt-RPL36-S19 Rabbit pAb
Species	Human
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S19 of alt-RPL36
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human
Shipping	Ice bag