

Phospho-Alt-Rpl36-S22 Rabbit Polyclonal Antibody - Alt-Rpl36

Item Number: CM0001167

Introduction

Phospho-alt-RPL36-S22 rabbit pAb for WB, DB, ELISA. Detects phosphorylated alt-RPL36 at Ser22 in human samples. Ideal for studying ribosomal protein phosphorylation dynamics.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-alt-RPL36-S22 Rabbit pAb
Species	Human
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around S22 of alt-RPL36
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human
Shipping	Ice bag