

Phospho-Ampk α 1-S485 Rabbit Polyclonal Antibody For Wb, Elisa - Q13131

Item Number: CM0001311



Introduction

Phospho-AMPK α 1-S485 rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, rat AMPK α 1 phosphorylated at S485, ~64kDa. Ideal for AMPK signaling and energy metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AMPK α 1-S485 Rabbit pAb
Remarks/Alias	AMPK; AMPK α 1; AMPK alpha 1; Phospho-AMPK α 1-S485
Species	Human
Gene ID (Human)	5562
Gene ID	5562
Immunogen	A synthetic phosphorylated peptide around S485 of human PRKAA1 (NP_006242.5).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13131
Protein Weight	64kDa
Shipping	Ice bag