

Phospho-Ampk α 1-S496 Rabbit Monoclonal Antibody For Wb, Elisa - Q13131

Item Number: CM0001645



Introduction

Phospho-specific rabbit monoclonal antibody recognizing AMPK α 1 phosphorylated at Ser496. Suitable for Western blot and ELISA applications. Validated in human and mouse samples. Ships on ice.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AMPK α 1-S496 Rabbit mAb
Remarks/Alias	AMPK; AMPK α 1; AMPK alpha 1; Phospho-AMPK α 1-S496
Species	Human
Gene ID (Human)	5562
Gene ID	5562
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S496 of human AMPK α 1 (Q13131).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q13131
Protein Weight	64kDa
Shipping	Ice bag