

Phospho-AMPK α 1-T183+AMPK α 2-T172 Rabbit Mab - Q13131 / P54646

Item Number: CM0001562



Introduction

Phospho-AMPK α 1-T183+AMPK α 2-T172 Rabbit mAb for WB and ELISA. Detects activated AMPK α in human, mouse, rat samples. Supports research on energy metabolism, autophagy, and cancer signaling pathways.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-AMPK α 1-T183+AMPK α 2-T172 Rabbit mAb
Remarks/Alias	AMPK α 1/AMPK α 2; Phospho-AMPK α 1-T183+AMPK α 2-T172
Species	Human
Gene ID (Human)	5562/5563
Immunogen	Synthetic phosphorylated peptide around T183/T172 of human PRKAA1/PRKAA2 (NP_006242.5/NP_006243.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13131/P54646
Protein Weight	64kDa/65kDa/62kDa
Shipping	Ice bag