

Phospho-Ampka1-T183/Ampka2-T172 Rabbit Polyclonal Antibody For Wb, Elisa - Q13131 / P54646

Item Number: CM0000538



Introduction

Phospho-specific rabbit polyclonal antibody detecting AMPK α 1-T183 and AMPK α 2-T172 phosphorylation. Validated for WB and ELISA in human, mouse, and rat samples. Ideal for energy metabolism and AMPK signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AMPK α 1-T183/AMPK α 2-T172 Rabbit pAb
Remarks/Alias	AMPK α 1/AMPK α 2; Phospho-AMPK α 1-T183/AMPK α 2-T172
Species	Human
GeneID (Human)	5562/5563
GeneID	5562/5563
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T183 of human PRKAA1 (NP_006242.5).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13131/P54646
Protein Weight	64kDa/65kDa/62kDa
Shipping	Ice bag