

Phospho-AP2M1-T156 Rabbit Pab For Wb And Elisa - Q96Cw1

Item Number: CM0000520



Introduction

Phospho-specific rabbit polyclonal antibody against AP2M1 phosphorylated at Thr-156. Suitable for Western blotting and ELISA with human and mouse reactivity. Ideal for investigating clathrin-mediated endocytosis and AP-2 complex regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AP2M1-T156 Rabbit pAb
Remarks/Alias	mu2; AP50; MRD60; CLAPM1; Phospho-AP2M1-T156
Species	Human
Gene ID (Human)	1173
Gene ID	1173
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around T156 of human AP2M1 (NP_001020376.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96CW1
Protein Weight	50kDa
Shipping	Ice bag