

Phospho-APP-T668 Rabbit Pab For Wb And Elisa - P05067

Item Number: CM0001965



Introduction

Rabbit polyclonal antibody against phosphorylated APP (Thr668) for Western blot and ELISA. Detects endogenous phospho-APP in human, mouse, and rat samples. Ideal for Alzheimer's disease research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-APP-T668 Rabbit pAb
Remarks/Alias	AAA; AD1; PN2; ABPP; APPI; CVAP; ABETA; PN-II; preA4; CTFgamma; alpha-sAPP; Phospho-APP-T668
Species	Human
Gene ID (Human)	351
Gene ID	351
Immunogen	A synthetic phosphorylated peptide around T668 of human APP (NP_958817.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P05067
Protein Weight	87kDa
Shipping	Ice bag