


Phospho-Arhgef2-S886 Rabbit Polyclonal Antibody For Wb, Elisa - Q92974

Item Number: CM0001846



Introduction

Rabbit polyclonal antibody against phosphorylated ARHGEF2 at Serine 886, validated for WB and ELISA in human, mouse, and rat. Suitable for studying Rho GTPase signaling, cell migration, and innate immune response.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ARHGEF2-S886 Rabbit pAb
Remarks/Alias	GEF; Lfc; P40; GEFH1; LFP40; GEF-H1; NEDMHM; Phospho-ARHGEF2-S886
Species	Human
GeneID (Human)	9181
GeneID	9181
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S886 of human ARHGEF2 (NP_001155855.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92974
Protein Weight	112kDa
Shipping	Ice bag