

Phospho-Ask1-T845 Rabbit Mab For Wb, Db, Elisa - O35099

Item Number: CM0001929

Introduction

Phospho-specific rabbit monoclonal antibody detecting ASK1 phosphorylated at Thr845. Validated for WB, DB, ELISA in mouse; detects the 155 kDa stress-activated kinase involved in apoptosis and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ASK1-T845 Rabbit mAb
Remarks/Alias	ASK; ASK1; Mekk5; MAPKKK5; 7420452D20Rik
Species	Mouse
Gene ID (Human)	4217
Gene ID	26408
Immunogen	A synthetic phosphorylated peptide around T845 of mouse ASK1 (NP_032606.4).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Mouse
SWISS	O35099
Protein Weight	155kDa
Shipping	Ice bag