

Phospho-Atf-2 (Thr69/71)/Atf-7 (Thr51/53) Rabbit Polyclonal Antibody For Wb, Elisa - P15336/P17544

Item Number: CM0001230



Introduction

Rabbit polyclonal antibody recognizing phosphorylated ATF-2 (Thr69/71)/ATF-7 (Thr51/53). Validated for WB and ELISA in human samples, suitable for studying DNA damage response and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ATF-2 (Thr69/71)/ATF-7 (Thr51/53) Rabbit pAb
Remarks/Alias	HB16; CREB2; TREB7; CREB-2; CRE-BP1; ATFA; Phospho-ATF-2 (Thr69/71)/ATF-7 (Thr51/53)
Species	Human
GeneID (Human)	1386/11016
GeneID	1386/11016
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around Thr69/71 and Thr51/53 of human Phospho-ATF-2/ATF-7.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P15336/P17544
Protein Weight	54kDa/51kDa
Shipping	Ice bag