


Phospho-Atf2 (T69/T71) Rabbit Pab For Wb, Ihc-P, Elisa - P15336

Item Number: CM0000866



Introduction

Phospho-specific rabbit polyclonal antibody targeting human ATF2 at Thr69/Thr71. Validated for Western blot, IHC-P, and ELISA in human, mouse, rat. Ideal for studying apoptosis, transcription, and DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ATF2-T69/T71 Rabbit pAb
Remarks/Alias	HB16; CREB2; TREB7; CREB-2; CRE-BP1; Phospho-ATF2-T69/T71
Species	Human
GeneID (Human)	1386
GeneID	1386
Immunogen	A synthetic phosphorylated peptide around T69 & T71 of human ATF2 (NP_001871.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15336
Protein Weight	55kDa
Shipping	Ice bag