

Phospho-Atf2-T71 Rabbit Mab For Wb And Elisa - P15336

Item Number: CM0000690



Introduction

Phospho-ATF2-T71 Rabbit mAb detects ATF2 phosphorylated at Thr71, validated for WB and ELISA with mouse cross-reactivity. Ideal for transcription regulation and DNA damage studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ATF2-T71 Rabbit mAb
Remarks/Alias	HB16; CREB2; TREB7; CREB-2; CRE-BP1; Phospho-ATF2-T71
Species	Human
GeneID (Human)	1386
GeneID	1386
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P15336
Protein Weight	55kDa
Shipping	Ice bag