


Phospho-Atf2-T71 Rabbit Pab For Wb, Ihc-P, Elisa - P15336

Item Number: CM0000579



Introduction

High-quality rabbit polyclonal antibody specific for phosphorylated ATF2 at Thr71. Validated for WB, IHC-P, ELISA; cross-reacts with human, mouse, rat. Used for transcription and DNA damage research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ATF2-T71 Rabbit pAb
Remarks/Alias	HB16; CREB2; TREB7; CREB-2; CRE-BP1; Phospho-ATF2-T71
Species	Human
Gene ID (Human)	1386
Gene ID	1386
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around T71 of human ATF2 (NP_001871.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15336
Protein Weight	55kDa
Shipping	Ice bag