

Phospho-Atm-S1981 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q13315

Item Number: CM0001078



Introduction

Phospho-S1981 ATM rabbit monoclonal antibody, validated for WB, IF-P, ELISA applications. Detects human ATM phosphorylated at serine 1981. For DNA damage response and checkpoint signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ATM-S1981 Rabbit mAb
Remarks/Alias	AT1; ATA; ATC; ATD; ATE; ATDC; TEL1; TEL01; Phospho-ATM-S1981
Species	Human
GeneID (Human)	472
GeneID	472
Immunogen	A synthetic phosphorylated peptide around S1981 of human ATM (Q13315)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	Q13315
Protein Weight	351kDa
Shipping	Ice bag