


Phospho-Atr-S428 Rabbit Mab - Q13535

Item Number: CM0001558



Introduction

Phospho-specific rabbit monoclonal antibody targeting human ATR at S428, validated for WB, IF-P, IHC-P, ELISA. Essential for DNA damage checkpoint and nuclear envelope integrity research, with 301kDa protein weight.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ATR-S428 Rabbit mAb
Remarks/Alias	FRP1; MEC1; SCKL; FCTCS; SCKL1; Phospho-ATR-S428
Species	Human
GeneID (Human)	545
GeneID	545
Immunogen	Synthetic peptide, A synthetic phosphorylated peptide around S428 of human ATR (NP_001175.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q13535
Protein Weight	301kDa
Shipping	Ice bag