


Phospho-Atr-S428 Rabbit Monoclonal Antibody For Wb And Elisa - Q13535

Item Number: CM0001807



Introduction

Phospho-ATR-S428 Rabbit mAb validated for WB, ELISA. Detects human ATR phosphorylated at serine 428, a critical DNA damage checkpoint kinase. Ideal for genotoxic stress and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ATR-S428 Rabbit mAb
Remarks/Alias	FRP1; MEC1; SCKL; FCTCS; SCKL1; Phospho-ATR-S428
Species	Human
Gene ID	545
Immunogen	Synthetic peptide / A synthetic phosphorylated peptide around S428 of human ATR (NP_001175.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13535
Protein Weight	301 kDa
Shipping	Ice bag