


Phospho-Atr-S428 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - Q13535

Item Number: CM0001761



Introduction

Rabbit polyclonal antibody specific for ATR phosphorylated at Ser428. Validated in WB, IHC-P, IF/ICC, and ELISA with reactivity across human, mouse, and rat. An essential tool for DNA damage checkpoint and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ATR-S428 Rabbit pAb
Remarks/Alias	FRP1; MEC1; SCKL; FCTCS; SCKL1; Phospho-ATR-S428
Species	Human
GeneID (Human)	545
GeneID	545
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around S428 of human ATR (NP_001175.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13535
Protein Weight	301kDa
Shipping	Ice bag