


Phospho-Atr-T1989 Rabbit Pab For Wb, Ihc-P, Elisa - Q13535

Item Number: CM0000712



Introduction

Phospho-ATR-T1989 Rabbit pAb detects phosphorylated ATR at T1989. Validated for WB, IHC-P, ELISA in human, mouse, rat. Essential for DNA damage response, checkpoint signaling, and genome stability research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ATR-T1989 Rabbit pAb
Remarks/Alias	FRP1; MEC1; SCKL; FCTCS; SCKL1; Phospho-ATR-T1989
Species	Human
GeneID (human)	545
GeneID	545
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13535
Protein Weight	301kDa
Shipping	Ice bag