


Phospho-Aurka (Thr288) Rabbit Polyclonal Antibody For Wb And Elisa - O14965

Item Number: CM0000889



Introduction

Rabbit polyclonal antibody specific for phosphorylated Aurora kinase A at Thr288. Validated for WB and ELISA in human samples. Suitable for mitosis, centrosome function, and cancer signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AURKA-T288 Rabbit pAb
Remarks/Alias	AIK; ARK1; AURA; BTAK; STK6; STK7; STK15; PPP1R47; Phospho-AURKA-T288
Species	Human
GeneID (human)	6790
GeneID	6790
Immunogen	A synthetic phosphorylated peptide around T288 of human AURKA (NP_003591.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14965
Protein Weight	46kDa
Shipping	Ice bag