


Phospho-Axl (Y702) Rabbit Polyclonal Antibody For Wb, Elisa - P30530

Item Number: CM0001486



Introduction

Phospho-specific rabbit polyclonal antibody targeting Axl at tyrosine 702 (P30530). Validated in WB and ELISA, suitable for studying GAS6/AXL signaling, cell survival, and immune response.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-Axl-Y702 Rabbit pAb
Remarks/Alias	ARK; UFO; JTK11; Tyro7; Phospho-Axl-Y702
Species	Human
Gene ID (Human)	558
Gene ID	558
Immunogen	Synthetic phosphorylated peptide around Y702 of human AXL (NP_068713.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P30530
Protein Weight	98kDa
Shipping	Ice bag