

Phospho-B-Raf (S445) Polyclonal Antibody For Wb, Elisa - P15056

Item Number: CM0001775



Introduction

Phospho-B-Raf (S445) Rabbit pAb detects endogenous B-Raf only when phosphorylated at S445. Validated for WB and ELISA with human and rat reactivity. Suitable for MAP kinase pathway studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-B-Raf-S445 Rabbit pAb
Remarks/Alias	NS7; B-raf; BRAF1; RAFB1; B-RAF1; BRAF-1; Phospho-B-Raf-S445
Species	Human
GeneID (Human)	673
GeneID	673
Immunogen	Synthetic peptide; A synthetic phosphorylated peptide around S445 of human BRAF (NP_004324.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P15056
Protein Weight	84kDa
Shipping	Ice bag