

Phospho-B-Raf-S728 Rabbit Pab For Wb - P15056

Item Number: CM0001593



Introduction

Phospho-B-Raf-S728 Rabbit pAb is a rabbit polyclonal antibody targeting human B-Raf phosphorylated at Ser728. Validated for WB and ELISA, it detects an 84kDa band. Ideal for studying MAPK signaling and oncogenic BRAF mutations.

[Learn More](#)

Parameter	Value
Product Name	Phospho-B-Raf-S728 Rabbit pAb
Remarks/Alias	NS7; B-raf; BRAF1; RAFB1; B-RAF1; BRAF-1; Phospho-B-Raf-S728
Species	Human
GeneID (human)	673
GeneID	673
Immunogen	A synthetic phosphorylated peptide around S728 of human BRAF (NP_004324.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P15056
Protein Weight	84kDa
Shipping	Ice bag