

Phospho-Bad-S112 Rabbit Mab For Wb, Ihc-P, Elisa - Q92934

Item Number: CM0001827



Introduction

Phospho-Bad-S112 Rabbit mAb is a rabbit monoclonal antibody targeting S112-phosphorylated human Bad, validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat Bad. Ideal for apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Bad-S112 Rabbit mAb
Remarks/Alias	BBC2; BCL2L8; Phospho-Bad-S112
Species	Human
Gene ID	572
Immunogen	Synthetic Peptide A phospho specific peptide corresponding to residues surrounding S112 of human Bad.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92934
Protein Weight	18kDa
Shipping	Ice bag