

Phospho-Bcar1-Y165 Rabbit Pab For Wb, Elisa - P56945

Item Number: CM0000982



Introduction

Phospho-specific BCAR1 (pY165) rabbit polyclonal antibody validated for Western blot and ELISA. Detects human, mouse, and rat p130cas. Ideal for studying integrin-mediated signaling, cell adhesion, and migration.

[Learn More](#)

Parameter	Value
Product Name	Phospho-BCAR1-Y165 Rabbit pAb
Remarks/Alias	CAS, CAS1, CASS1, CRKAS, P130Cas, Phospho-BCAR1-Y165
Species	Human
Gene ID (Human)	9564
Gene ID	9564
Immunogen	Synthetic Peptide - A synthetic phosphorylated peptide around Y165 of human BCAR1 (NP_055382.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56945
Protein Weight	93kDa
Shipping	Ice bag