

Phospho-Bcar1 (Y410) Rabbit Polyclonal Antibody For Wb, Elisa - P56945

Item Number: CM0001452



Introduction

Rabbit polyclonal antibody targeting phosphorylated BCAR1 at tyrosine 410 (P56945). Validated for WB and ELISA, reacts with human and mouse samples. Ideal for cell adhesion and migration signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-BCAR1-Y410 Rabbit pAb
Remarks/Alias	CAS; CAS1; CASS1; CRKAS; P130Cas; Phospho-BCAR1-Y410
Species	Human
GeneID (Human)	9564
GeneID	9564
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y410 of human BCAR1 (NP_055382.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P56945
Protein Weight	93kDa
Shipping	Ice bag