

Phospho-Bcl-2-S87 Rabbit Pab For Wb - P10415

Item Number: CM0001094



Introduction

Phospho-specific rabbit polyclonal antibody detecting Bcl-2 phosphorylated at Ser87. Ideal for WB and ELISA in human and rat samples; supports apoptosis, autophagy, and cancer research.

[Learn More](#)

Parameter	Details
Product Name	Phospho-Bcl-2-S87 Rabbit pAb
Remarks/Alias	Bcl-2; PPP1R50; Phospho-Bcl-2-S87
Species	Human
GeneID (Human)	596
GeneID	596
Immunogen	A synthetic phosphorylated peptide around S87 of human Bcl-2 (NP_000624.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P10415
Protein Weight	26kDa
Shipping	Ice bag