

## Phospho-Beclin 1 (Ser15) Rabbit Polyclonal Antibody For Wb And Elisa - Q14457

Item Number: CM0001642



### Introduction

Phospho-specific Beclin-1 (Ser15) rabbit polyclonal antibody for WB and ELISA. Reacts with human, mouse, and rat Beclin-1 (52 kDa). Essential for autophagy, apoptosis, and host defense research.

[Learn More](#)

| Parameter        | Specification                                                                  |
|------------------|--------------------------------------------------------------------------------|
| Product Name     | Phospho-Beclin 1-Ser15 Rabbit pAb                                              |
| Remarks/Alias    | ATG6; VPS30; beclin1; Phospho-Beclin 1-Ser15                                   |
| Species          | Human                                                                          |
| Gene ID (Human)  | 8678                                                                           |
| Gene ID          | 8678                                                                           |
| Immunogen        | A synthetic phosphorylated peptide around S15 of human Beclin 1 (NP_003757.1). |
| Source           | Rabbit                                                                         |
| Category         | Phosphorylated Antibodies                                                      |
| Application      | WB, ELISA                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                              |
| SWISS            | Q14457                                                                         |
| Protein Weight   | 52kDa                                                                          |
| Shipping         | Ice bag                                                                        |