

Phospho-Beclin 1 (Ser90) Rabbit Polyclonal Antibody For Wb And Elisa - Q14457

Item Number: CM0001314



Introduction

Phospho-specific rabbit polyclonal antibody targeting Beclin-1 at Ser90. Validated for Western blot and ELISA in human, mouse, and rat. Suitable for autophagy, apoptosis, and endocytosis research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Beclin 1-Ser90 Rabbit pAb
Remarks/Alias	ATG6; VPS30; beclin1; Phospho-Beclin 1-Ser90
Species	Human
GeneID (Human)	8678
GeneID	8678
Immunogen	A synthetic phosphorylated peptide around S90 of human Beclin 1 (NP_003757.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14457
Protein Weight	52kDa
Shipping	Ice bag