

Phospho-Beclin 1-Ser93 Rabbit Pab For Wb, Elisa - Q14457

Item Number: CM0001425



Introduction

Phospho-Beclin 1-Ser93 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human Beclin-1 phosphorylated at Ser93, a key autophagy and apoptosis regulator. Ideal for research in autophagy, apoptosis, and host defense.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Beclin 1-Ser93 Rabbit pAb
Remarks/Alias	ATG6; VPS30; beclin1; Phospho-Beclin 1-Ser93
Species	Human
Gene ID (Human)	8678
Gene ID	8678
Immunogen	A synthetic phosphorylated peptide around S93 of human Beclin 1 (NP_003757.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14457
Protein Weight	52kDa
Shipping	Ice bag