

Phospho-Brca1-S1423 Rabbit Pab - P38398

Item Number: CM0000740



Introduction

Phospho-BRCA1 (S1423) rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA with human, mouse, rat reactivity. Useful for DNA damage repair and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-BRCA1-S1423 Rabbit pAb
Remarks/Alias	IRIS; PSCP; BRCAI; BRCC1; FANCS; PNCA4; RNF53; BROVCA1; PPP1R53; Phospho-BRCA1-S1423
Species	Human
GeneID (Human)	672
GeneID	672
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S1423 of human BRCA1 (NP_009225.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P38398
Protein Weight	208kDa
Shipping	Ice bag