


Phospho-Brd4 (S492/S494) Rabbit Polyclonal Antibody For Wb, Db, Elisa - O60885

Item Number: CM0000726



Introduction

Phospho-specific rabbit polyclonal antibody targeting BRD4 at S492/S494. Validated for Western blot, dot blot, and ELISA with human reactivity; studies epigenetic memory and transcription regulation. ~152 kDa. UniProt O60885.

[Learn More](#)

Parameter	Value
Product Name	Phospho-BRD4-S492/S494 Rabbit pAb
Remarks/Alias	CAP; MCAP; HUNK1; HUNKI; FSHRG4; Phospho-BRD4-S492/S494
Species	Human
GeneID (human)	23476
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human
SWISS	O60885
Protein Weight	152kDa
Shipping	Ice bag