


Phospho-Brd4-S498/S499 Rabbit Polyclonal Antibody For Wb, Elisa - O60885

Item Number: CM0000845



Introduction

This rabbit polyclonal antibody detects endogenous levels of BRD4 only when phosphorylated at S498/S499. Suitable for WB and ELISA with human reactivity. Ideal for epigenetic memory and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-BRD4-S498/S499 Rabbit pAb
Remarks/Alias	CAP; MCAP; HUNK1; HUNKI; FSHRG4; Phospho-BRD4-S498/S499
Species	Human
GeneID	23476
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60885
Protein Weight	152kDa
Shipping	Ice bag