


## Phospho-Brd4-T204 Rabbit Mab For Wb, Elisa - O60885

Item Number: CM0001163



### Introduction

Phospho-BRD4-T204 Rabbit monoclonal antibody validated for WB and ELISA. Detects human and rat BRD4 phosphorylated at Threonine 204. Suitable for studying epigenetic memory, transcription regulation, and NF-kappa-B signaling.

[Learn More](#)

| Parameter        | Value                                                                      |
|------------------|----------------------------------------------------------------------------|
| Product Name     | Phospho-BRD4-T204 Rabbit mAb                                               |
| Remarks/Alias    | CAP; MCAP; HUNK1; HUNKI; FSHRG4; Phospho-BRD4-T204                         |
| Species          | Human                                                                      |
| GeneID (Human)   | 23476                                                                      |
| GeneID           | 23476                                                                      |
| Immunogen        | A synthetic phosphorylated peptide around T204 of human BRD4 (NP_490597.1) |
| Source           | Rabbit                                                                     |
| Category         | Phosphorylated Antibodies                                                  |
| Application      | WB, ELISA                                                                  |
| Cross Reactivity | Human,Rat                                                                  |
| SWISS            | O60885                                                                     |
| Protein Weight   | 152kDa                                                                     |
| Shipping         | Ice bag                                                                    |