


Phospho-Brd4-T204 Rabbit Pab - O60885

Item Number: CM0000823



Introduction

Phospho-BRD4-T204 Rabbit pAb detects phosphorylated BRD4 at threonine 204. Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Suitable for epigenetic and transcription regulation studies. ~152kDa.

[Learn More](#)

Parameter	Value
Product Name	Phospho-BRD4-T204 Rabbit pAb
Remarks/Alias	CAP; MCAP; HUNK1; HUNKI; FSHRG4; Phospho-BRD4-T204
Species	Human
GeneID (Human)	23476
GeneID	23476
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T204 of human BRD4 (NP_490597.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60885
Protein Weight	152kDa
Shipping	Ice bag