


Phospho-Brd4-T249 Rabbit Pab For Wb - O60885

Item Number: CM0001199



Introduction

Phospho-BRD4-T249 Rabbit pAb detects BRD4 phosphorylated at Thr249. Validated for WB in human, mouse, and rat samples. Ideal for epigenetic memory and transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-BRD4-T249 Rabbit pAb
Remarks/Alias	CAP; MCAP; HUNK1; HUNKI; FSHRG4; Phospho-BRD4-T249
Species	Human
Gene ID (Human)	23476
Gene ID	23476
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around T249 of human BRD4 (NP_490597.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human, Mouse, Rat
SWISS	O60885
Protein Weight	152kDa
Shipping	Ice bag