


Phospho-C-Abl-Y412 Rabbit Pab - P00519

Item Number: CM0000491



Introduction

Phospho-specific rabbit polyclonal antibody detecting c-Abl phosphorylated at Y412. Suitable for IHC-P and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for studying ABL1 kinase signaling, apoptosis, and cytoskeletal dynamics.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Abl-Y412 Rabbit pAb
Remarks/Alias	ABL; JTK7; p150; c-ABL; v-abl; CHDSKM; c-ABL1; BCR-ABL; bcr/abl; Phospho-c-Abl-Y412
Species	Human
GeneID (human)	25
GeneID	25
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y412 of human ABL1 (NP_009297.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00519
Protein Weight	123kDa
Shipping	Ice bag