


Phospho-C-Abl-Y412 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P00519

Item Number: CM0000426



Introduction

Rabbit polyclonal antibody specific for phospho-c-Abl (Y412) for use in WB, IHC-P, and ELISA. Detects c-Abl phosphorylation in human, mouse, and rat samples. Ideal for cancer, kinase signaling, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Abl-Y412 Rabbit pAb
Remarks/Alias	ABL; JTK7; p150; c-ABL; v-abl; CHDSKM; c-ABL1; BCR-ABL; bcr/abl; Phospho-c-Abl-Y412
Species	Human
GeneID (Human)	25
GeneID	25
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around Y412 of human ABL1 (NP_005148.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00519
Protein Weight	123kDa
Shipping	Ice bag