


Phospho-C-Jun-S243 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P05412

Item Number: CM0000599



Introduction

Phospho-c-Jun-S243 Rabbit pAb detects human, mouse, and rat c-Jun phosphorylated at Ser243. Validated for WB, IHC-P, and ELISA, this antibody supports AP-1 transcription complex and cancer signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Jun-S243 Rabbit pAb
Remarks/Alias	AP1; p39; AP-1; cJUN; c-Jun; Phospho-c-Jun-S243
Species	Human
Gene ID (Human)	3725
Gene ID	3725
Immunogen	Synthetic peptide / A synthetic phosphorylated peptide around S243 of human Jun (NP_002219.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05412
Protein Weight	36kDa
Shipping	Ice bag