


Phospho-C-Jun-S63 Rabbit Pab For Wb And Elisa - P05412

Item Number: CM0000877



Introduction

Phospho-specific rabbit polyclonal antibody targeting c-Jun phosphorylated at Ser63. Validated for Western blot and ELISA, cross-reacts with human and rat. Ideal for studying AP-1 transcription factor signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Jun-S63 Rabbit pAb
Remarks/Alias	AP1; p39; AP-1; cJUN; c-Jun; Phospho-c-Jun-S63
Species	Human
GeneID (Human)	3725
GeneID	3725
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S63 of human c-Jun (NP_002219.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P05412
Protein Weight	36kDa
Shipping	Ice bag