


Phospho-C-Jun-S63 Rabbit Pab - P05412

Item Number: CM0001709



Introduction

Phospho-c-Jun-S63 Rabbit pAb detects phosphorylated c-Jun at Ser63 in human, mouse, and rat. Suitable for WB, IHC-P, ELISA in transcription and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Jun-S63 Rabbit pAb
Remarks/Alias	AP1; p39; AP-1; cJUN; c-Jun; Phospho-c-Jun-S63
Species	Human
GeneID (Human)	3725
GeneID	3725
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S63 of human Jun (NP_002219.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05412
Protein Weight	36kDa
Shipping	Ice bag