


Phospho-C-Jun-S73 Rabbit Pab - P05412

Item Number: CM0001691



Introduction

Rabbit polyclonal antibody targeting phosphorylated c-Jun at Serine 73, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Ideal for studying AP-1 transcription factor signaling in cancer and immune responses.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Jun-S73 Rabbit pAb
Remarks/Alias	AP1; p39; AP-1; cJUN; c-Jun; Phospho-c-Jun-S73
Species	Human
GeneID (Human)	3725
GeneID	3725
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S73 of human c-Jun (NP_002219.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05412
Protein Weight	36kDa
Shipping	Ice bag