


Phospho-C-Jun-S73 Rabbit Pab - P05412

Item Number: CM0001064



Introduction

Phospho-c-Jun-S73 Rabbit polyclonal antibody for WB, IHC-P, ELISA. Detects human, mouse, rat c-Jun phosphorylated at S73. IVD raw material for cancer and transcription research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Jun-S73 Rabbit pAb
Remarks/Alias	AP1; p39; AP-1; cJUN; c-Jun; Phospho-c-Jun-S73
Species	Human
GeneID	3725
Immunogen	A synthetic phosphorylated peptide around S73 of human c-JUN (NP_002219.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05412
Protein Weight	36kDa
Shipping	Ice bag