


Phospho-C-Jun-T91 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P05412

Item Number: CM0001220



Introduction

Phospho-specific rabbit polyclonal antibody targeting Thr91 of transcription factor c-Jun (AP-1, p39). Validated for WB, IHC-P, IF/ICC, and ELISA; reacts with human, mouse, and rat. Ideal for studying JNK signaling, cell proliferation, and oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Jun-T91 Rabbit pAb
Remarks/Alias	AP1; p39; AP-1; cJUN; c-Jun; Phospho-c-Jun-T91
Species	Human
Gene ID	3725
Immunogen	A synthetic phosphorylated peptide around T91 of human JUN (NP_002219.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05412
Protein Weight	36kDa
Shipping	Ice bag