


Phospho-C-Myc-S62 Rabbit Monoclonal Antibody For Wb And Elisa - P01106

Item Number: CM0001143



Introduction

High-specificity Phospho-c-Myc-S62 rabbit monoclonal antibody validated for Western blotting (WB) and ELISA in human and mouse samples. Ideal for studying c-Myc phosphorylation and transcriptional activity. Molecular weight ~51 kDa.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Myc-S62 Rabbit mAb
Remarks/Alias	MRTL; MYCC; c-Myc; bHLHe39; Phospho-c-Myc-S62
Species	Human
GeneID (Human)	4609
GeneID	4609
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S62 of human c-Myc (NP_002458.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01106
Protein Weight	51kDa
Shipping	Ice bag