


Phospho-C-Myc-T58/S62 Rabbit Mab For Wb, Elisa - P01106

Item Number: CM0000874



Introduction

Phospho-c-Myc (T58/S62) Rabbit Monoclonal Antibody validated for WB and ELISA. Specific to human c-Myc phosphorylated at Thr58 and Ser62, with rat cross-reactivity. Ideal for studying Myc signaling, cell proliferation, and oncogene activation.

[Learn More](#)

Specification	Details
Product Name	Phospho-c-Myc-T58/S62 Rabbit mAb
Remarks/Alias	MRTL; MYCC; c-Myc; bHLHe39; Phospho-c-Myc-T58/S62
Species	Human
Gene ID	4609
Immunogen	A synthetic phosphorylated peptide around T58 & S62 of human c-Myc (P01106)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P01106
Protein Weight	51kDa
Shipping	Ice bag