


Phospho-Cardiac Troponin I (Tnni3)-S23/S24 Rabbit Mab For Ihc-P, Elisa - P19429

Item Number: CM0001792



Introduction

Phospho-specific rabbit monoclonal antibody targeting phosphorylated S23/S24 of human cardiac troponin I (TNNI3), validated for IHC-P and ELISA. Detects human cardiac troponin I, a key regulator of cardiac muscle contraction. Ideal for cardiac biomarker and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Cardiac Troponin I (TNNI3)-S23/S24 Rabbit mAb
Remarks/Alias	CMH7; RCM1; cTnI; CMD2A; TNNC1; CMD1FF; Phospho-Cardiac Troponin I (TNNI3)-S23/S24
Species	Human
Gene ID (Human)	7137
Gene ID	7137
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S23 & S24 of human Cardiac Troponin I (TNNI3) (NP_000354.4).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P19429
Protein Weight	24kDa
Shipping	Ice bag