

Phospho-Casein Kinase 2 Alpha (Csnk2A1)-T360/S362 Rabbit Polyclonal Antibody For Wb And Elisa - P68400

Item Number: CM0000881



Introduction

Rabbit polyclonal antibody specific for CSNK2A1 phosphorylated at T360/S362. Validated for WB and ELISA, detects human and mouse Casein Kinase 2 alpha. Ideal for studying kinase signaling, apoptosis, and transcription.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Casein Kinase 2 alpha (CSNK2A1)-T360/S362 Rabbit pAb
Remarks/Alias	CKII; Cka1; Cka2; CK2A1; OCNDS; Phospho-Casein Kinase 2 alpha (CSNK2A1)-T360/S362
Species	Human
Gene ID (Human)	1457
Immunogen	Synthetic peptide: A phospho specific peptide corresponding to residues surrounding T360/S362 of human Casein Kinase 2 alpha (CSNK2A1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P68400
Protein Weight	45kDa
Shipping	Ice bag