


Phospho-Cbl-Y700 Rabbit Polyclonal Antibody For Wb And Elisa - P22681

Item Number: CM0000461



Introduction

Phospho-CBL-Y700 Rabbit Polyclonal Antibody detects phosphorylated CBL at Tyr700 in human samples, validated for WB and ELISA. Suitable for investigating ubiquitination, RTK signaling, and osteoclast biology.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CBL-Y700 Rabbit pAb
Remarks/Alias	CBL2; NSLL; C-CBL; RNF55; FRA11B; Phospho-CBL-Y700
Species	Human
Gene ID (Human)	867
Gene ID	867
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y700 of human CBL (NP_005179.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P22681
Protein Weight	100kDa
Shipping	Ice bag